


ProPrep Disc Preparation Set. For Anterior Lumbar Surgery.

Comprehensive set
of instruments

Facilitates
discectomy for
anterior lumbar
procedures

Includes reduced-
profile, angled,
and longer length
instruments




ProPrep Disc Preparation Set for Anterior Lumbar Surgery

A premium set of core instruments optimized for anterior lumbar interbody fusion and nonfusion procedures using an anterior approach.

Features

- Consists of a simplified but comprehensive set of instruments
- Facilitates a complete discectomy for anterior lumbar fusion and nonfusion procedures
- Addresses collapsed segments with reduced-profile instruments
- Allows access to posterolateral corners using angled instruments
- Provides optimized length for anterior approach and large body masses
- Incorporates “two-hand-grip” silicone handles for maximum instrument control



ProPrep Instruments

-
- 03.605.000 Disc Rongeurs, 330 mm length
Up-biting, (side teeth), 3 mm width
 - 03.605.001 Straight, without teeth, 4 mm width
 - 03.605.002 Straight, without teeth, 6 mm width



03.605.000

03.605.001

-
- 03.605.003 Box Curette, 6 mm x 10 mm,
430 mm length



-
- 03.605.004 Endplate Elevator, 20 mm width,
430 mm length



-
- Bone Curettes
angled, oval head, 430 mm length
 - 03.605.006 3.5 mm x 4.5 mm
 - 03.605.007 5.5 mm x 8.5 mm
 - 03.605.008 7.5 mm x 11.5 mm



03.605.006



03.605.007



03.605.008

03.605.010 Ball Tip Probe, 90°, 300 mm length



03.605.011 Ring Curette, teardrop, (angled),
8 mm width, 430 mm length



Laminectomy Punches
40° up-biting, 330 mm length
389.42 2 mm width
389.41 4 mm width
389.412 6 mm width



389.481 Bone Rongeur, double-acting without
teeth, 8 mm, 360 mm length




Also Available

PDL114 Vertebral Body Spreader, angled



03.605.005 Bone Elevator, 17 mm tip width



03.605.009 Midline Marker, 8 mm width



03.605.012 Blunt Dissector, 267 mm length



ProPrep Disc Preparation Set for Anterior Lumbar Surgery (01.605.600)

Graphic Case

690.210 Graphic Case for ProPrep Instrument Set

Instruments

Disc Rongeurs, 330 mm length

03.605.000 Up-biting, (side teeth), 3 mm width

03.605.001 Straight, without teeth, 4 mm width

03.605.002 Straight, without teeth, 6 mm width

03.605.003 Box Curette, 6 mm x 10 mm

03.605.004 Endplate Elevator, 20 mm width

Bone Curettes, angled, oval head

03.605.006 3.5 mm x 4.5 mm

03.605.007 5.5 mm x 8.5 mm

03.605.008 7.5 mm x 11.5 mm

03.605.010 Ball Tip Probe, 90°, 300 mm length

03.605.011 Ring Curette, teardrop, (angled), 8 mm width

40° Up-biting Laminectomy Punches, 330 mm length

389.42 2 mm width

389.41 4 mm width

389.412 6 mm width

389.481 Bone Rongeur, double-acting without teeth, 8 mm width, 360 mm length



Optional Instruments

PDL114 Vertebral Body Spreader, angled

03.605.005 Bone Elevator, 17 mm tip width

03.605.009 Midline Marker, 8 mm width, 250 mm length

03.605.012 Blunt Dissector, 267 mm

The Graphic Case (690.210) is designed to hold additional instruments and to allow substitution of similar instruments within the set, for procedures using an anterior approach.

Note: For additional information, please refer to package insert.

Also Available

	90° Up-biting Laminectomy Punches, 330 mm length
389.43	4 mm width
389.44	2 mm width
	Disc Rongeurs, up-biting, with teeth
389.46	4 mm width
389.47	7.5 mm width
389.56	Endplate Elevator, 25 mm width
389.57	Intervertebral Ring Curette, angled, 11 mm width
	Bone Curettes, angled, oval head
389.59	9.5 mm x 15 mm
389.591	11.5 mm x 17.5 mm
389.60	13.5 mm x 20 mm
389.61	3.5 mm x 4.5 mm
	Bone Curettes, straight, oval head
389.65	9.5 mm x 15 mm
389.651	11.5 mm x 17.5 mm
389.66	13.5 mm x 20 mm
389.67	3.5 mm x 4.5 mm
389.671	5.5 mm x 8 mm
389.672	7.5 mm x 11.5 mm
	Osteotomes, straight
389.73	18 mm width
389.74	25 mm width
389.75	Impactor, 12 mm x 20 mm



Synthes Spine
1302 Wrights Lane East
West Chester, PA 19380
Telephone: (610) 719-5000
To order: (800) 523-0322
Fax: (610) 251-9056

Synthes (Canada) Ltd.
2566 Meadowpine Boulevard
Mississauga, Ontario L5N 6P9
Telephone: (905) 567-0440
To order: (800) 668-1119
Fax: (905) 567-3185

www.synthes.com