

Large Oscillating Saw Attachment and Blades.

For Tibial Plateau Leveling Osteotomy (TPLO).

One power system for the entire surgery

Specially designed saw attachment for TPLO surgery

Saw blades for TPLO ranging from 12 mm to 30 mm



Tibial Plateau Leveling Osteotomy (TPLO)

Small Battery Drive System

- Compact battery driven power tool
- Light and ergonomic
- One power system for most orthopaedic applications, including drilling, screwing, sawing and K-wire insertion
- 14.4 V battery is recommended for use with the Oscillating Saw Attachment, large (532.026)
- Compact case size fits most tabletop autoclaves



Saw Attachment and Saw Blades for TPLO

- Specially designed saw attachment for tibial plateau leveling osteotomies and dome osteotomies
- Seven saw blades, with radii ranging from 12 mm to 30 mm
- Simple and stable connection of the saw blade into the saw attachment
- Thin saw blades offer excellent cutting performance and minimal bone removal (cutting thickness is 0.6 mm)
- Compatible with Synthes TPLO plate system



Product Information

Small Battery Drive and Accessories

530.601	Universal Battery Charger, for Synthes Batteries, 4 bays
532.004	Battery Insertion Shield
532.010	Small Battery Drive
532.026	Oscillating Saw Attachment, large
532.027	Replacement Key, for 532.026
532.032	Battery Casing for 14.4 V Battery
532.033	14.4 V Battery for Small Battery Drive



Crescent Saw Blades

03.000.390	12 mm radius, 45 mm x 0.6 mm
03.000.391	15 mm radius, 45 mm x 0.6 mm
03.000.392	18 mm radius, 45 mm x 0.6 mm
03.000.393*	21 mm radius, 45 mm x 0.6 mm
03.000.394	24 mm radius, 45 mm x 0.6 mm
03.000.395	27 mm radius, 50 mm x 0.6 mm
03.000.396	30 mm radius, 50 mm x 0.6 mm



03.000.390



03.000.391



03.000.392



03.000.394



03.000.395



03.000.396

All sawblades come with attachment screws.

* Not shown

For the full range of attachments and accessories for the small battery drive, please contact your Synthes representative or consult the Synthes product catalog.



Synthes (USA)
1302 Wrights Lane East
West Chester, PA 19380
Telephone: (610) 719-5000
To order: (800) 523-0322
Fax: (610) 251-9056

Synthes (Canada) Ltd.
2566 Meadowpine Boulevard
Mississauga, Ontario L5N 6P9
Telephone: (905) 567-0440
To order: (800) 668-1119
Fax: (905) 567-3185

www.synthes.com